

Product Datasheet

NLRP3/NALP3 Antibody [DyLight 650] AF7010W

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/AF7010W

Updated 5/28/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/AF7010W



AF7010W**NLRP3/NALP3 Antibody [DyLight 650]**

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 650
Purity	Antigen Affinity-purified
Buffer	50mM Sodium Borate
Product Description	
Host	Goat
Gene ID	114548
Gene Symbol	NLRP3
Species	Human
Specificity/Sensitivity	Detects human NLRP3/NALP3 in direct ELISAs. In direct ELISAs, less than 1% cross-reactivity with recombinant human (rh) NALP1, rhNALP2, and rhNALP5 is observed.
Immunogen	<i>E. coli</i> -derived recombinant human NLRP3/NALP3 Lys540-Lys689 Accession # Q96P20
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Immunocytochemistry
Recommended Dilutions	Immunocytochemistry
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to AF7010W

NB410-28088C-500ul	Goat IgG Isotype Control [DyLight 650]
H00114548-P01-2ug	Recombinant Human NLRP3/NALP3 GST (N-Term) Protein
201-LB-005	IL-1 beta/IL-1F2 [Unconjugated]
NBP1-77080PEP	NLRP3/NALP3 Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/AF7010W

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

