

# Product Datasheet

## Katanin p60 Antibody [Biotin] AF7100B

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF7100B](http://www.novusbio.com/AF7100B)

Updated 5/28/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF7100B](http://www.novusbio.com/reviews/destination/AF7100B)



**AF7100B**

Katanin p60 Antibody [Biotin]

| Product Information         |                                                                                              |
|-----------------------------|----------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                       |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.  |
| Storage                     | Store at 4C in the dark.                                                                     |
| Clonality                   | Polyclonal                                                                                   |
| Preservative                | 0.05% Sodium Azide                                                                           |
| Isotype                     | IgG                                                                                          |
| Conjugate                   | Biotin                                                                                       |
| Purity                      | Antigen Affinity-purified                                                                    |
| Buffer                      | PBS                                                                                          |
| Product Description         |                                                                                              |
| Host                        | Sheep                                                                                        |
| Gene ID                     | 11104                                                                                        |
| Gene Symbol                 | KATNA1                                                                                       |
| Species                     | Human                                                                                        |
| Specificity/Sensitivity     | Detects human Katanin p60 in direct ELISAs.                                                  |
| Immunogen                   | <i>E. coli</i> -derived recombinant human Katanin p60<br>Asp308-Tyr378<br>Accession # O75449 |
| Product Application Details |                                                                                              |
| Applications                | Immunohistochemistry                                                                         |
| Recommended Dilutions       | Immunohistochemistry                                                                         |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                       |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF7100B**

---

|               |                                           |
|---------------|-------------------------------------------|
| BAF020        | Sheep IgG Isotype Control [Biotin]        |
| NBP1-89491PEP | Katanin p60 Recombinant Protein Antigen   |
| NBP2-13075    | beta Amyloid Antibody (MOAB-2) - BSA Free |
| NBP2-01383    | katanin-p80 Antibody (OTI5C4)             |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF7100B](http://www.novusbio.com/reviews/submit/AF7100B)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

