

# Product Datasheet

## Spi-B Antibody [PerCP] AF7204C

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF7204C](http://www.novusbio.com/AF7204C)

Updated 5/28/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF7204C](http://www.novusbio.com/reviews/destination/AF7204C)



**AF7204C**

Spi-B Antibody [PerCP]

| Product Information         |                                                                                                                                 |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                          |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                     |
| Storage                     | Store at 4C in the dark.                                                                                                        |
| Clonality                   | Polyclonal                                                                                                                      |
| Preservative                | 0.05% Sodium Azide                                                                                                              |
| Isotype                     | IgG                                                                                                                             |
| Conjugate                   | PerCP                                                                                                                           |
| Purity                      | Antigen Affinity-purified                                                                                                       |
| Buffer                      | PBS                                                                                                                             |
| Product Description         |                                                                                                                                 |
| Host                        | Sheep                                                                                                                           |
| Gene ID                     | 6689                                                                                                                            |
| Gene Symbol                 | SPIB                                                                                                                            |
| Species                     | Mouse                                                                                                                           |
| Specificity/Sensitivity     | Detects mouse Spi-B in direct ELISAs. In direct ELISAs, less than 2% cross-reactivity with recombinant human Spi-B is observed. |
| Immunogen                   | <i>E. coli</i> -derived recombinant mouse Spi-B<br>Tyr18-Glu167 (Tyr110Phe)<br>Accession # O35906                               |
| Product Application Details |                                                                                                                                 |
| Applications                | CyTOF-ready, Immunocytochemistry, Intracellular Staining by Flow Cytometry                                                      |
| Recommended Dilutions       | Intracellular Staining by Flow Cytometry, Immunocytochemistry, CyTOF-ready                                                      |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                          |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF7204C**

---

|                    |                                              |
|--------------------|----------------------------------------------|
| H00006689-P01-10ug | Recombinant Human Spi-B GST (N-Term) Protein |
| 210-TA-005         | TNF-alpha [Unconjugated]                     |
| NBL1-16399         | Spi-B Overexpression Lysate                  |
| DY417-05           | IL-10 [Biotin]                               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF7204C](http://www.novusbio.com/reviews/submit/AF7204C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

