

# Product Datasheet

## Lymphotoxin-alpha/TNF-beta Antibody [DyLight 488] AF749K

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF749K](http://www.novusbio.com/AF749K)

Updated 7/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF749K](http://www.novusbio.com/reviews/destination/AF749K)



**AF749K****Lymphotoxin-alpha/TNF-beta Antibody [DyLight 488]**

| <b>Product Information</b>         |                                                                                                                                                                                                   |
|------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 ml                                                                                                                                                                                            |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                       |
| <b>Storage</b>                     | Store at 4C in the dark.                                                                                                                                                                          |
| <b>Clonality</b>                   | Polyclonal                                                                                                                                                                                        |
| <b>Preservative</b>                | 0.05% Sodium Azide                                                                                                                                                                                |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                               |
| <b>Conjugate</b>                   | DyLight 488                                                                                                                                                                                       |
| <b>Purity</b>                      | Antigen Affinity-purified                                                                                                                                                                         |
| <b>Buffer</b>                      | 50mM Sodium Borate                                                                                                                                                                                |
| <b>Product Description</b>         |                                                                                                                                                                                                   |
| <b>Host</b>                        | Goat                                                                                                                                                                                              |
| <b>Gene ID</b>                     | 4049                                                                                                                                                                                              |
| <b>Gene Symbol</b>                 | LTA                                                                                                                                                                                               |
| <b>Species</b>                     | Mouse                                                                                                                                                                                             |
| <b>Specificity/Sensitivity</b>     | Detects mouse Lymphotoxin $\alpha$ /TNF $\beta$ in direct ELISAs and Western blots. In direct ELISAs, 100% cross-reactivity with recombinant human Lymphotoxin $\alpha$ /TNF $\beta$ is observed. |
| <b>Immunogen</b>                   | <i>E. coli</i> -derived recombinant mouse Lymphotoxin $\alpha$ /TNF $\beta$                                                                                                                       |
| <b>Notes</b>                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                 |
| <b>Product Application Details</b> |                                                                                                                                                                                                   |
| <b>Applications</b>                | Western Blot                                                                                                                                                                                      |
| <b>Recommended Dilutions</b>       | Western Blot                                                                                                                                                                                      |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                            |



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF749K**

|                 |                                                      |
|-----------------|------------------------------------------------------|
| NB410-28088G    | Goat IgG Isotype Control [DyLight 488]               |
| NBP2-61307-10ug | Recombinant Human Lymphotoxin-alpha/TNF-beta Protein |
| 210-TA-005      | TNF-alpha [Unconjugated]                             |
| 211-TBB-010/CF  | Lymphotoxin-alpha/TNF-beta [Unconjugated]            |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF749K](http://www.novusbio.com/reviews/submit/AF749K)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

