

# Product Datasheet

## Caspase-9 Antibody [PE] AF8301P

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF8301P](http://www.novusbio.com/AF8301P)

Updated 5/8/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF8301P](http://www.novusbio.com/reviews/destination/AF8301P)



**AF8301P****Caspase-9 Antibody [PE]**

| <b>Product Information</b>         |                                                                                                                                                                                                                                                                                                                                   |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 ml                                                                                                                                                                                                                                                                                                                            |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                       |
| <b>Storage</b>                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                          |
| <b>Clonality</b>                   | Polyclonal                                                                                                                                                                                                                                                                                                                        |
| <b>Preservative</b>                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                               |
| <b>Conjugate</b>                   | PE                                                                                                                                                                                                                                                                                                                                |
| <b>Purity</b>                      | Antigen Affinity-purified                                                                                                                                                                                                                                                                                                         |
| <b>Buffer</b>                      | PBS                                                                                                                                                                                                                                                                                                                               |
| <b>Product Description</b>         |                                                                                                                                                                                                                                                                                                                                   |
| <b>Host</b>                        | Goat                                                                                                                                                                                                                                                                                                                              |
| <b>Gene ID</b>                     | 842                                                                                                                                                                                                                                                                                                                               |
| <b>Gene Symbol</b>                 | CASP9                                                                                                                                                                                                                                                                                                                             |
| <b>Species</b>                     | Human                                                                                                                                                                                                                                                                                                                             |
| <b>Specificity/Sensitivity</b>     | Detects human Caspase-9 precursor and the pro + large subunit (LSU) that migrates as a 37 kDa band on SDS-PAGE generated by cleavage at Asp330. Does not detect the pro + LSU that migrates as a 34 kDa band on SDS-PAGE generated by cleavage at Asp315. Therefore, the major epitope detected occurs between Asp315 and Asp330. |
| <b>Immunogen</b>                   | <i>E. coli</i> -derived recombinant human Caspase-9<br>Val139-Asp330<br>Accession # P55211                                                                                                                                                                                                                                        |
| <b>Product Application Details</b> |                                                                                                                                                                                                                                                                                                                                   |
| <b>Applications</b>                | Western Blot                                                                                                                                                                                                                                                                                                                      |
| <b>Recommended Dilutions</b>       | Western Blot                                                                                                                                                                                                                                                                                                                      |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                            |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to AF8301P**

---

|                     |                                          |
|---------------------|------------------------------------------|
| NB410-28088PE-500ul | Goat IgG Isotype Control [PE]            |
| NB100-56121PEP      | Caspase-9 Antibody Blocking Peptide      |
| 210-TA-005          | TNF-alpha [Unconjugated]                 |
| NBP2-75042          | Mouse Caspase-9 ELISA Kit (Colorimetric) |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF8301P](http://www.novusbio.com/reviews/submit/AF8301P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

