

Product Datasheet

CD40/TNFRSF5 Antibody (mAb89) - Azide and BSA Free DDX0131P-100

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 13

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/DDX0131P-100

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/DDX0131P-100



DDX0131P-100


CD40/TNFRSF5 Antibody (mAb89) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	mAb89
Preservative	No Preservative
Isotype	IgG1
Purity	Ion exchange chromatography
Buffer	Tris-NaCl (pH 8.0)
Product Description	
Host	Mouse
Gene ID	958
Gene Symbol	CD40
Species	Human
Immunogen	Anti-u activated tonsillar human B cells.
Notes	This product is manufactured by Eurobio Scientific (formerly Dendritics) and distributed by Novus Biologicals.
Product Application Details	
Applications	Flow Cytometry, Functional, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin
Recommended Dilutions	Flow Cytometry 1:10-1:1000, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence 1:10-1:500, Immunohistochemistry-Paraffin 1:10-1:500, Immunohistochemistry-Frozen 1:10-1:500, Functional




Images

Immunocytochemistry/Immunofluorescence: CD40/TNFRSF5 Antibody (mAb89) [DDX0131P-100] - CD40-transfected L cells stained.



Immunohistochemistry: CD40/TNFRSF5 Antibody (mAb89) [DDX0131P-100] - IHC staining on Bouin paraffin lymph node section. Human CD40 stained with mAb 89.



Publications

Warner BD, Carbognani A, Franco L et Al. LIGHTCURVE ANALYSIS OF THE NEAR-EARTH ASTEROID 2015 TB145 Minor Planet Bull 2016-01-01 [PMID: 32455372]

Warner BD, Harris AW, Pravec P et Al. LIGHTCURVE PHOTOMETRY OPPORTUNITIES: 2014 JULY-SEPTEMBER Minor Planet Bull 2014-01-01 [PMID: 32457962]

de Saint Vis B et al. A novel lysosome-associated membrane glycoprotein, DC-LAMP, induced upon DC maturation, is transiently expressed in MHC class II compartment. Immunity 1998-09-01 [PMID: 9768752]

Bjorck P et al. CD40 antibodies defining distinct epitopes display qualitative differences in their induction of B-cell differentiation. Immunology 1996-02-01 [PMID: 8698393]

Bjorck P et al. Antibodies to distinct epitopes on the CD40 molecule co-operate in stimulation and can be used for the detection of soluble CD40. Immunology 1994-11-01 [PMID: 7530692]

Lozza L et al. Crosstalk between human DC subsets promotes antibacterial activity and CD8+ T-cell stimulation in response to bacille Calmette-Guerin. Eur J Immunol 2014-01-01 [PMID: 24114554]

Fromm JR et al. Flow cytometric analysis of CD123 is useful for immunophenotyping classical Hodgkin lymphoma. Cytometry B Clin Cytom 2011-03-01 [PMID: 20809502]

Meyer-Wentrup F et al. Targeting DCIR on human plasmacytoid dendritic cells results in antigen presentation and inhibits IFN-alpha production. Blood 2008-04-01 [PMID: 18258799]

Contin C et al. Membrane-anchored CD40 is processed by the tumor necrosis factor-alpha-converting enzyme. Implications for CD40 signaling. J Biol Chem 2003-08-01 [PMID: 12810728]

Valladeau J et al. Immature human dendritic cells express asialoglycoprotein receptor isoforms for efficient receptor-mediated endocytosis. J Immunol 2001-11-01 [PMID: 11698450]

Henn V et al. The inflammatory action of CD40 ligand (CD154) expressed on activated human platelets is temporally limited by coexpressed CD40. Blood 2001-08-01 [PMID: 11493450]

Tong AW et al. Growth-inhibitory effects of CD40 ligand (CD154) and its endogenous expression in human breast cancer. Clin Cancer Res 2001-03-01 [PMID: 11297266]

More publications at <http://www.novusbio.com/DDX0131P-100>





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to DDX0131P-100

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP2-51628-0.05mg	Recombinant Human CD40/TNFRSF5 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/DDX0131P-100

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

