

# Product Datasheet

## Human IgG3 Isotype Control [Alexa Fluor® 488] DDXCH03A488-100

Unit Size: 0.1 mg

Store at -20C in the dark. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/DDXCH03A488-100](http://www.novusbio.com/DDXCH03A488-100)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/DDXCH03A488-100](http://www.novusbio.com/reviews/destination/DDXCH03A488-100)



**DDXCH03A488-100**

Human IgG3 Isotype Control [Alexa Fluor® 488]

<b>Product Information</b>	
<b>Unit Size</b>	0.1 mg
<b>Concentration</b>	Please contact technical services for concentration.
<b>Storage</b>	Store at -20C in the dark. Avoid freeze-thaw cycles.
<b>Clonality</b>	Polyclonal
<b>Preservative</b>	No Preservative
<b>Isotype</b>	IgG3
<b>Conjugate</b>	Alexa Fluor 488
<b>Purity</b>	Ion exchange chromatography
<b>Buffer</b>	PBS, 50% glycerol.

<b>Product Description</b>	
<b>Description</b>	This Human IgG3 immunoglobulin is useful as an isotype control for Human IgG3 antibody.
<b>Host</b>	Human
<b>Gene ID</b>	3502
<b>Gene Symbol</b>	IGHG3
<b>Species</b>	Human
<b>Notes</b>	This product is manufactured by Eurobio Scientific (formerly Dendritics) and distributed by Novus Biologicals.

<b>Product Application Details</b>	
<b>Applications</b>	Flow Cytometry, Flow (Cell Surface), Flow (Intracellular), Immunofluorescence
<b>Recommended Dilutions</b>	Flow Cytometry 1:10-1:1000, Flow (Cell Surface) 1:10 - 1:1000, Flow (Intracellular) 1:10 - 1:1000, Immunofluorescence 1:10 - 1:500
<b>Application Notes</b>	Use the isotype control at the same concentration as the experimental antibody. Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to DDXCH03A488-100**

---

H00003502-P01-10ug	Recombinant Human IgG3 GST (N-Term) Protein
210-TA-005	TNF-alpha [Unconjugated]
NBP2-62125	IgG3 Antibody (RM119)
M6000B-1	IL-6 [HRP]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Isotype Controls are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/DDXCH03A488-100](http://www.novusbio.com/reviews/submit/DDXCH03A488-100)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

