

Product Datasheet

Mouse IgG3 Isotype Control [Alexa Fluor® 647] DDXCM03A647

Unit Size: 0.1 mg

Store at -20C in the dark. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/DDXCM03A647

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/DDXCM03A647



DDXCM03A647

Mouse IgG3 Isotype Control [Alexa Fluor® 647]

Product Information	
Unit Size	0.1 mg
Concentration	Please contact technical services for concentration.
Storage	Store at -20C in the dark. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG3
Conjugate	Alexa Fluor 647
Purity	Ion exchange chromatography
Buffer	PBS, 50% glycerol.

Product Description	
Description	This Mouse IgG3 immunoglobulin is useful as an isotype control for Mouse IgG3 antibody.
Host	Mouse
Gene ID	3502
Gene Symbol	IGHG3
Species	Mouse
Notes	This product is manufactured by Eurobio Scientific (formerly Dendritics) and distributed by Novus Biologicals.

Product Application Details	
Applications	Flow (Cell Surface), Flow (Intracellular), Immunofluorescence
Recommended Dilutions	Flow (Cell Surface) 1:10 - 1:1000, Flow (Intracellular) 1:10 - 1:1000, Immunofluorescence 1:10 - 1:500
Application Notes	Use the isotype control at the same concentration as the experimental antibody. Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to DDXCM03A647

H00003502-P01-10ug	Recombinant Human IgG3 GST (N-Term) Protein
210-TA-005	TNF-alpha [Unconjugated]
NBP2-62125	IgG3 Antibody (RM119)
M6000B-1	IL-6 [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Isotype Controls are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/DDXCM03A647

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

