

Product Datasheet

p23/PTGES3 Antibody (998320) [DyLight 488] FAB10039K

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB10039K

Updated 5/8/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB10039K



FAB10039K

p23/PTGES3 Antibody (998320) [DyLight 488]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	998320
Preservative	0.05% Sodium Azide
Isotype	IgG2a
Conjugate	DyLight 488
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	50mM Sodium Borate

Product Description	
Host	Mouse
Gene ID	10728
Species	Human, Mouse, Rat
Specificity/Sensitivity	Detects human p23/PTGES3 in direct ELISAs. Detects human, mouse, and rat p23/PTGES3 in Western blots.
Immunogen	<i>E. coli</i> -derived recombinant Human p23/PTGES3 Gln21-Glu160 Accession # Q15185
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB10039K

NBP3-18323-100ug	Recombinant Human p23/PTGES3 Protein
NBL1-14937	p23/PTGES3 Overexpression Lysate
NB200-106	p19ARF/CDKN2A Antibody - BSA Free
NB100-689	COX-2 Antibody

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB10039K

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

