

Product Datasheet

Lymphotoxin beta R/TNFRSF3 Antibody (157108)

[PE/Cy5.5]

FAB10081PECY55

Unit Size: 0.5 ml

Store at 4C in the dark. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB10081PECY55

Updated 8/9/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB10081PECY55



FAB10081PECY55

Lymphotoxin beta R/TNFRSF3 Antibody (157108) [PE/Cy5.5]

| Product Information | |
|-----------------------------|--|
| Unit Size | 0.5 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. Do not freeze. |
| Clonality | Monoclonal |
| Clone | 157108 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2a |
| Conjugate | PE/Cy5.5 |
| Purity | Protein A or G purified from hybridoma culture supernatant |
| Buffer | PBS |
| Product Description | |
| Host | Rat |
| Gene ID | 4055 |
| Gene Symbol | LTBR |
| Species | Mouse |
| Specificity/Sensitivity | Detects mouse Lymphotoxin β R in direct ELISAs and Western blots. In direct ELISAs, does not cross-react with recombinant mouse (rm) 4 \square 1BB, rmBAFF R, rmCD27, rmCD30, rmCD40, recombinant human (rh) DR3, rhDR6, rmEDAR, rmFas, rmGITR, rmHVEM, rhLymphotoxin β R, rmNGF R, rmOPG, rmOX40, rmRANK, rhRELT, rmTAJ, rmTNF RI, rmTNF RII, rhTRADD, rmTRAIL R2, rhTRAIL R3, rhTRAIL R4, or rmTWEAK R. |
| Immunogen | Mouse myeloma cell line NS0-derived recombinant mouse Lymphotoxin β R/TNFRSF3 Ser28-Pro218 Accession # P50284 |
| Product Application Details | |
| Applications | Flow Cytometry |
| Recommended Dilutions | Flow Cytometry |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. For optimal results using our Tandem dyes, please avoid prolonged exposure to light or extreme temperature fluctuations. These can lead to irreversible degradation or decoupling. When staining intracellular targets, specific attention to the fixation and permeabilization steps in your flow protocol may be required. Please contact our technical support team at technical@novusbio.com if you have any questions. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB10081PECY55

| | |
|-------------------|--|
| NBP2-51906-0.05mg | Recombinant Human Lymphotoxin beta R/TNFRSF3 His Protein |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| 7538-LR-100 | Lymphotoxin beta R/TNFRSF3 [Unconjugated] |
| M6000B-1 | IL-6 [HRP] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB10081PECY55

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

