

# Product Datasheet

## Betacellulin/BTC Antibody (197415) [DyLight 680] FAB1025Y

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB1025Y](http://www.novusbio.com/FAB1025Y)

Updated 7/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB1025Y](http://www.novusbio.com/reviews/destination/FAB1025Y)



**FAB1025Y**

Betacellulin/BTC Antibody (197415) [DyLight 680]

| Product Information         |                                                                                                                                                                                                                                                                                                                 |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.5 ml                                                                                                                                                                                                                                                                                                          |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                     |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                        |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                      |
| Clone                       | 197415                                                                                                                                                                                                                                                                                                          |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                              |
| Isotype                     | IgG2a                                                                                                                                                                                                                                                                                                           |
| Conjugate                   | DyLight 680                                                                                                                                                                                                                                                                                                     |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                                      |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                              |
| Product Description         |                                                                                                                                                                                                                                                                                                                 |
| Host                        | Rat                                                                                                                                                                                                                                                                                                             |
| Gene ID                     | 685                                                                                                                                                                                                                                                                                                             |
| Gene Symbol                 | BTC                                                                                                                                                                                                                                                                                                             |
| Species                     | Mouse                                                                                                                                                                                                                                                                                                           |
| Specificity/Sensitivity     | Detects Betacellulin/BTC in direct ELISAs and Western blots. In direct ELISAs and Western blots, no cross-reactivity with recombinant human (rh) Androgen R, rhBetacellulin, rhEGF, rhHB-EGF, rhHRG- $\alpha$ , rhHRG- $\beta$ 1, rhTGF- $\alpha$ , recombinant mouse (rm) Epigen, or rmEpiregulin is observed. |
| Immunogen                   | <i>E. coli</i> -derived recombinant mouse Betacellulin/BTC<br>Asp32-Gln118<br>Accession # Q05928                                                                                                                                                                                                                |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                               |
| Product Application Details |                                                                                                                                                                                                                                                                                                                 |
| Applications                | Western Blot                                                                                                                                                                                                                                                                                                    |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                                                                                                                                    |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                          |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB1025Y**

---

|                |                                            |
|----------------|--------------------------------------------|
| NBP2-35148-5ug | Recombinant Mouse Betacellulin/BTC Protein |
| 236-EG-200     | EGF [Unconjugated]                         |
| 261-CE-010     | Betacellulin/BTC [Unconjugated]            |
| 338-AC-010     | Activin A [Unconjugated]                   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB1025Y](http://www.novusbio.com/reviews/submit/FAB1025Y)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

