

# Product Datasheet

## Kallikrein 5 Antibody (325302) [DyLight 405] FAB11082E

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB11082E](http://www.novusbio.com/FAB11082E)

Updated 5/28/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB11082E](http://www.novusbio.com/reviews/destination/FAB11082E)



**FAB11082E**

Kallikrein 5 Antibody (325302) [DyLight 405]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                         |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                    |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                       |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                                     |
| Clone                       | 325302                                                                                                                                                                                                                                                                                                                                                         |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                             |
| Isotype                     | IgG1                                                                                                                                                                                                                                                                                                                                                           |
| Conjugate                   | DyLight 405                                                                                                                                                                                                                                                                                                                                                    |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                                                                                     |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                                             |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                                          |
| Gene ID                     | 25818                                                                                                                                                                                                                                                                                                                                                          |
| Gene Symbol                 | KLK5                                                                                                                                                                                                                                                                                                                                                           |
| Species                     | Human                                                                                                                                                                                                                                                                                                                                                          |
| Specificity/Sensitivity     | Detects human Kallikrein 5 Propeptide in direct ELISAs and Western blots. In direct ELISAs and Western blots, no cross-reactivity with recombinant human (rh) Kallikrein 5 (mature), rhKallikrein 1, -3, -4, -6, -7, -8, -9, -10, -11, -12, -13, -14, -15, rhKallikrein B1, rhHGFA, rhFactor VII, rhFactor X, rhFactor XI, rhThrombin, or rhuPA-1 is observed. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant human Kallikrein 5 Val23-Arg66<br>Accession # NP_036559                                                                                                                                                                                                                                                        |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                                              |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                |
| Applications                | Western Blot                                                                                                                                                                                                                                                                                                                                                   |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                                                                                                                                                                                   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                         |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB11082E**

---

|                   |                                                |
|-------------------|------------------------------------------------|
| NBP1-97005V-0.5ml | Mouse IgG1 Isotype Control (MG1) [DyLight 405] |
| NBP2-59548-50ug   | Recombinant Human Kallikrein 5 His Protein     |
| 1398-TS-010       | TSLP [Unconjugated]                            |
| 1108-SE-010       | Kallikrein 5 [Unconjugated]                    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB11082E](http://www.novusbio.com/reviews/submit/FAB11082E)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

