

# Product Datasheet

## Dkk-3 Antibody (138728) [FITC] FAB11181F

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB11181F](http://www.novusbio.com/FAB11181F)

Updated 7/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB11181F](http://www.novusbio.com/reviews/destination/FAB11181F)



**FAB11181F**

Dkk-3 Antibody (138728) [FITC]

| Product Information         |                                                                                                                                                          |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.5 ml                                                                                                                                                   |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                              |
| Storage                     | Store at 4C in the dark.                                                                                                                                 |
| Clonality                   | Monoclonal                                                                                                                                               |
| Clone                       | 138728                                                                                                                                                   |
| Preservative                | 0.05% Sodium Azide                                                                                                                                       |
| Isotype                     | IgG2a                                                                                                                                                    |
| Conjugate                   | FITC                                                                                                                                                     |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                               |
| Buffer                      | PBS                                                                                                                                                      |
| Product Description         |                                                                                                                                                          |
| Host                        | Mouse                                                                                                                                                    |
| Gene ID                     | 27122                                                                                                                                                    |
| Gene Symbol                 | DKK3                                                                                                                                                     |
| Species                     | Human                                                                                                                                                    |
| Specificity/Sensitivity     | Detects human Dkk-3 in ELISAs. In direct ELISAs, no cross-reactivity with recombinant human Dkk-1, recombinant mouse (rm) Dkk-2, or rmDkk-3 is observed. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant human Dkk-3<br>Ala22-Ile350<br>Accession # Q9UBP4.1                                                      |
| Product Application Details |                                                                                                                                                          |
| Applications                | ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), ELISA Standard (Matched Pair)                                            |
| Recommended Dilutions       | ELISA Standard (Matched Pair), ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair)                                            |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                   |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB11181F**

---

|               |                                          |
|---------------|------------------------------------------|
| NBP1-96848    | Mouse IgG2a Isotype Control (M2A) [FITC] |
| NBP1-80843PEP | Dkk-3 Recombinant Protein Antigen        |
| 5439-DK-010   | Dkk-1 [Unconjugated]                     |
| 1118-DK-050   | Dkk-3 [Unconjugated]                     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB11181F](http://www.novusbio.com/reviews/submit/FAB11181F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

