

# Product Datasheet

## 5-HT2A Antibody (1055220) [DyLight 550] FAB112501L

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB112501L](http://www.novusbio.com/FAB112501L)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB112501L](http://www.novusbio.com/reviews/destination/FAB112501L)



**FAB112501L**

5-HT2A Antibody (1055220) [DyLight 550]

| Product Information |   |
|---------------------|---|
| Unit Size           | 0.1 ml  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Monoclonal  |
| Clone               | 1055220   |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG1  |
| Conjugate           | DyLight 550   |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | 50mM Sodium Borate  |

| Product Description     |  |
|-------------------------|--|
| Host                    | Mouse  |
| Gene ID                 | 3356   |
| Species                 | Human  |
| Specificity/Sensitivity | Detects human HTR2A in direct ELISA.   |
| Immunogen               | Synthetic peptide corresponding to the amino acids 1-75 in the N-term extracellular domain of the Serotonin Receptor 5-HT2A (HTR2A).<br>Accession # P28223 |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.  |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | Flow Cytometry, Immunocytochemistry                                    |
| Recommended Dilutions       | Flow Cytometry, Immunocytochemistry                                    |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB112501L**

---

|             |  |
|-------------|--|
| NBP1-97005R | Mouse IgG1 Isotype Control (MG1) [DyLight 550] |
| NBP3-31022  | Human 5-HT2A ELISA Kit (Colorimetric)          |
| NB100-56352 | 5-HT7 Antibody                                 |
| NB100-1533  | POMC Antibody                                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB112501L](http://www.novusbio.com/reviews/submit/FAB112501L)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

