

# Product Datasheet

## **E6AP/UBE3A Antibody (1077315) [DyLight 650] FAB11546W**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB11546W](http://www.novusbio.com/FAB11546W)

Updated 11/12/2025 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB11546W](http://www.novusbio.com/reviews/destination/FAB11546W)



**FAB11546W**

E6AP/UBE3A Antibody (1077315) [DyLight 650]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 1077315                                                                                     |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG1                                                                                        |
| Conjugate           | DyLight 650                                                                                 |
| Purity              | Protein A or G purified from cell culture supernatant                                       |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description     |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| Host                    | Mouse                                                                                             |
| Gene ID                 | 7337                                                                                              |
| Species                 | Human                                                                                             |
| Specificity/Sensitivity | Detects recombinant human His6-E6AP/UBE3A in Direct ELISA.                                        |
| Immunogen               | Spodoptera frugiperda Sf 21 (baculovirus)-derived human E6AP/UBE3A Met1-Leu875 Accession # Q05086 |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                 |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Immunohistochemistry, Immunocytochemistry                              |
| Recommended Dilutions       | Immunohistochemistry, Immunocytochemistry                              |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB11546W**

|                    |                                                   |
|--------------------|---------------------------------------------------|
| NBP1-97005C        | Mouse IgG1 Isotype Control (MG1) [DyLight 650]    |
| H00007337-P01-10ug | Recombinant Human E6AP/UBE3A GST (N-Term) Protein |
| E2-640-100         | Ubch7/UBE2L3 [Unconjugated]                       |
| E3-230-050         | E6AP/UBE3A [Unconjugated]                         |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB11546W](http://www.novusbio.com/reviews/submit/FAB11546W)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

