

# Product Datasheet

## IL-18 R beta/IL-1 R7/ACPL Antibody (132016) [PerCP] FAB1181C

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB1181C](http://www.novusbio.com/FAB1181C)

Updated 7/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB1181C](http://www.novusbio.com/reviews/destination/FAB1181C)



**FAB1181C**

IL-18 R beta/IL-1 R7/ACPL Antibody (132016) [PerCP]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.5 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 132016                                                                                      |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG1                                                                                        |
| Conjugate           | PerCP                                                                                       |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | PBS                                                                                         |

| Product Description     |                                                                                                                                                                                                                                                                                        |
|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Mouse                                                                                                                                                                                                                                                                                  |
| Gene ID                 | 8807                                                                                                                                                                                                                                                                                   |
| Gene Symbol             | IL18RAP                                                                                                                                                                                                                                                                                |
| Species                 | Human                                                                                                                                                                                                                                                                                  |
| Specificity/Sensitivity | Detects human IL-18 R $\beta$ /IL-1 R7 in direct ELISAs and Western blots. In direct ELISAs and Western blots, less than 1% cross-reactivity with recombinant mouse IL-18 R $\beta$ , recombinant human (rh) IL-18 R, rhIL-1 RI, rhIL-1 RII, rhIL-1 RAcP, and rhIL-1 Rrp2 is observed. |
| Immunogen               | Mouse myeloma cell line NS0-derived recombinant human IL-18 R $\beta$ /IL-1 R7 Met1-Arg356<br>Accession # O95256                                                                                                                                                                       |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot, Neutralization                                           |
| Recommended Dilutions       | Western Blot, Neutralization                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB1181C**

---

|                |                                                     |
|----------------|-----------------------------------------------------|
| NBP1-97005PCP  | Mouse IgG1 Isotype Control (MG1) [PerCP]            |
| NB100-56433PEP | IL-18 R beta/IL-1 R7/ACPL Antibody Blocking Peptide |
| 210-TA-005     | TNF-alpha [Unconjugated]                            |
| AVI118-050     | IL-18 R beta/IL-1 R7/ACPL [Biotin]                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB1181C](http://www.novusbio.com/reviews/submit/FAB1181C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

