

# Product Datasheet

## **Follistatin-related Gene Protein/FLRG/Fstl3 Antibody (206714) [DyLight 650] FAB1288W**

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB1288W](http://www.novusbio.com/FAB1288W)

Updated 7/14/2025 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB1288W](http://www.novusbio.com/reviews/destination/FAB1288W)



**FAB1288W**

Follistatin-related Gene Protein/FLRG/Fstl3 Antibody (206714) [DyLight 650]

| Product Information         |                                                                                                                                                                                           |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.5 ml                                                                                                                                                                                    |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                               |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                  |
| Clonality                   | Monoclonal                                                                                                                                                                                |
| Clone                       | 206714                                                                                                                                                                                    |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                        |
| Isotype                     | IgG1                                                                                                                                                                                      |
| Conjugate                   | DyLight 650                                                                                                                                                                               |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                        |
| Product Description         |                                                                                                                                                                                           |
| Host                        | Mouse                                                                                                                                                                                     |
| Gene ID                     | 10272                                                                                                                                                                                     |
| Gene Symbol                 | FSTL3                                                                                                                                                                                     |
| Species                     | Human                                                                                                                                                                                     |
| Specificity/Sensitivity     | Detects human Follistatin-related Gene Protein/FLRG in ELISAs and Western blots. In sandwich immunoassays and Western blots, no cross-reactivity with recombinant mouse FLRG is observed. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant human Follistatin-related Gene Protein/FLRG<br>Met27-Val263<br>Accession # O95633                                                         |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                         |
| Product Application Details |                                                                                                                                                                                           |
| Applications                | ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), ELISA Standard (Matched Pair)                                                                             |
| Recommended Dilutions       | ELISA Standard (Matched Pair), ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair)                                                                             |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                    |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB1288W**

---

|                |                                                                         |
|----------------|-------------------------------------------------------------------------|
| NBP1-97005C    | Mouse IgG1 Isotype Control (MG1) [DyLight 650]                          |
| NBP2-32034PEP  | Follistatin-related Gene Protein/FLRG/Fstl3 Recombinant Protein Antigen |
| 7754-BH-005/CF | TGF-beta 1 [Unconjugated]                                               |
| DFLRG0         | Follistatin-related Gene Protein/FLRG/Fstl3 [HRP]                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB1288W](http://www.novusbio.com/reviews/submit/FAB1288W)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

