

# Product Datasheet

## Angiopoietin-like Protein 3/ANGPTL3 Antibody (128610) [DyLight 680] FAB136Y

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB136Y](http://www.novusbio.com/FAB136Y)

Updated 7/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB136Y](http://www.novusbio.com/reviews/destination/FAB136Y)



**FAB136Y**

Angiopoietin-like Protein 3/ANGPTL3 Antibody (128610) [DyLight 680]

| Product Information         |                                                                                                                                                                                                                                                                                                   |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.5 ml                                                                                                                                                                                                                                                                                            |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                       |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                          |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                        |
| Clone                       | 128610                                                                                                                                                                                                                                                                                            |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                |
| Isotype                     | IgG2a                                                                                                                                                                                                                                                                                             |
| Conjugate                   | DyLight 680                                                                                                                                                                                                                                                                                       |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                        |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                |
| Product Description         |                                                                                                                                                                                                                                                                                                   |
| Host                        | Rat                                                                                                                                                                                                                                                                                               |
| Gene ID                     | 27329                                                                                                                                                                                                                                                                                             |
| Gene Symbol                 | ANGPTL3                                                                                                                                                                                                                                                                                           |
| Species                     | Mouse                                                                                                                                                                                                                                                                                             |
| Specificity/Sensitivity     | Detects mouse Angiopoietin-like 3 in direct ELISAs and Western blots. In direct ELISAs, this antibody does not cross-react with recombinant human (rh) Angiopoietin-1, rhAngiopoietin-2, rhAngiopoietin-4, rmAngiopoietin-3, and cornea-derived transcript 6 (CDT6), an Angiopoietin-like factor. |
| Immunogen                   | <i>S. frugiperda</i> insect ovarian cell line Sf 21-derived recombinant mouse Angiopoietin-like 3<br>Ser17-Thr455<br>Accession # Q9R182                                                                                                                                                           |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                 |
| Product Application Details |                                                                                                                                                                                                                                                                                                   |
| Applications                | Western Blot                                                                                                                                                                                                                                                                                      |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                                                                                                                      |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                            |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB136Y**

---

|            |                                           |
|------------|-------------------------------------------|
| DANG20     | Angiopoietin-2 [HRP]                      |
| DANL30     | Angiopoietin-like Protein 3/ANGPTL3 [HRP] |
| 923-AN-025 | Angiopoietin-1 [Unconjugated]             |
| MAB1417    | Insulin Antibody (182410) [Unconjugated]  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB136Y](http://www.novusbio.com/reviews/submit/FAB136Y)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

