

# Product Datasheet

## Fetuin A/AHSG Antibody (207440) [DyLight 650] FAB15632W

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB15632W](http://www.novusbio.com/FAB15632W)

Updated 7/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB15632W](http://www.novusbio.com/reviews/destination/FAB15632W)



**FAB15632W**

Fetuin A/AHSG Antibody (207440) [DyLight 650]

| Product Information         |                                                                                                                                                                                                                                                            |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.5 ml                                                                                                                                                                                                                                                     |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                   |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                 |
| Clone                       | 207440                                                                                                                                                                                                                                                     |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                         |
| Isotype                     | IgG2b                                                                                                                                                                                                                                                      |
| Conjugate                   | DyLight 650                                                                                                                                                                                                                                                |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                 |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                         |
| Product Description         |                                                                                                                                                                                                                                                            |
| Host                        | Rat                                                                                                                                                                                                                                                        |
| Gene ID                     | 197                                                                                                                                                                                                                                                        |
| Gene Symbol                 | AHSG                                                                                                                                                                                                                                                       |
| Species                     | Mouse                                                                                                                                                                                                                                                      |
| Specificity/Sensitivity     | Detects mouse Fetuin A/AHSG in ELISAs. In sandwich immunoassays, no cross-reactivity or interference was observed with recombinant human Fetuin A, rhTGF- $\beta$ 1, recombinant mouse (rm) BMP-4, rmMMP-2, recombinant rat MMP-2, or rmMMP-9 is observed. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant mouse Fetuin A/AHSG Ala19-Ile345<br>Accession # P29699                                                                                                                                                     |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                          |
| Product Application Details |                                                                                                                                                                                                                                                            |
| Applications                | ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), ELISA Standard (Matched Pair)                                                                                                                                              |
| Recommended Dilutions       | ELISA Standard (Matched Pair), ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair)                                                                                                                                              |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                     |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB15632W**

---

|               |                                           |
|---------------|-------------------------------------------|
| NBP1-90302PEP | Fetuin A/AHSG Recombinant Protein Antigen |
| 210-TA-005    | TNF-alpha [Unconjugated]                  |
| DFTA00        | Fetuin A/AHSG [HRP]                       |
| M6000B-1      | IL-6 [HRP]                                |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB15632W](http://www.novusbio.com/reviews/submit/FAB15632W)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

