

Product Datasheet

IL-12/IL-35 p35 Antibody (27537) [mFluor Violet 500 SE] FAB1570MFV500

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB1570MFV500

Updated 10/10/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB1570MFV500



FAB1570MFV500

IL-12/IL-35 p35 Antibody (27537) [mFluor Violet 500 SE]

| Product Information | |
|---------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 27537 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 |
| Conjugate | mFluor Violet 500 SE |
| Purity | 0 |
| Buffer | 50mM Sodium Borate |

| Product Description | |
|-------------------------|--|
| Host | Mouse |
| Gene ID | 3592 |
| Gene Symbol | IL12A |
| Species | Human, Mouse |
| Specificity/Sensitivity | Detects human IL \square 12/IL-35 p35 in direct ELISAs and Western blots. Detects the p35 subunit either as part of a p40/p35 heterodimer or as a free subunit after reduction of the heterodimer. This antibody does not recognize IL-12 p40 homodimers but shows strong cross-reactivity with the p35 subunits from porcine and mouse systems. |
| Immunogen | <i>S. frugiperda</i> insect ovarian cell line Sf 21-derived recombinant human IL \square 12/IL-35 p35 Arg23-Ser219 Accession # P29459 |
| Notes | mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |

| Product Application Details | |
|-----------------------------|--|
| Applications | Western Blot, CyTOF-ready, Intracellular Staining by Flow Cytometry |
| Recommended Dilutions | Western Blot, Intracellular Staining by Flow Cytometry, CyTOF-ready |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB1570MFV500

| | |
|------------|--|
| 210-TA-005 | TNF-alpha [Unconjugated] |
| NBP3-06742 | Human IL-12/IL-35 p35 ELISA Kit (Colorimetric) |
| 285-IF-100 | IFN-gamma [Unconjugated] |
| AF835 | Caspase-3 Antibody [Unconjugated] - Active |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB1570MFV500

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

