

# Product Datasheet

## TWEAK R/TNFRSF12 Antibody (314615) [HRP] FAB1610H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB1610H](http://www.novusbio.com/FAB1610H)

Updated 5/28/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB1610H](http://www.novusbio.com/reviews/destination/FAB1610H)



**FAB1610H**

TWEAK R/TNFRSF12 Antibody (314615) [HRP]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 314615                                                                                      |
| Preservative        | No Preservative                                                                             |
| Isotype             | IgG2b                                                                                       |
| Conjugate           | HRP                                                                                         |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | PBS                                                                                         |

| Product Description     |                                                                                                                                                                                                                                                                                                                                                               |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Rat                                                                                                                                                                                                                                                                                                                                                           |
| Gene ID                 | 51330                                                                                                                                                                                                                                                                                                                                                         |
| Gene Symbol             | TNFRSF12A                                                                                                                                                                                                                                                                                                                                                     |
| Species                 | Mouse                                                                                                                                                                                                                                                                                                                                                         |
| Specificity/Sensitivity | Detects mouse TWEAK R/TNFRSF12 in direct ELISAs and Western blots. In Western blots, no cross-reactivity with recombinant mouse (rm) 4-1BB, rmTNF RII, rmTRAIL R2, rmTNF RI, rmFAS, rmEDAR, rmRANK, rmCD30, rmOPG, rmDR3, rmNGF R, rmLTR $\beta$ , rmHVEM, rmGITR, rmCD27, rmBAFF, recombinant human (rh) DR6, rhRELt, rhTRAIL R3, or rhTRAIL R4 is observed. |
| Immunogen               | <i>E. coli</i> -derived recombinant mouse TWEAK R/TNFRSF12 Glu28-Trp79<br>Accession # Q9CR75                                                                                                                                                                                                                                                                  |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot                                                           |
| Recommended Dilutions       | Western Blot                                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB1610H**

---

|                |                                            |
|----------------|--------------------------------------------|
| NBP2-35085-5ug | Recombinant Human TWEAK R/TNFRSF12 Protein |
| 210-TA-005     | TNF-alpha [Unconjugated]                   |
| 1610-TW-050    | TWEAK R/TNFRSF12 [Unconjugated]            |
| DVE00          | VEGF [HRP]                                 |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB1610H](http://www.novusbio.com/reviews/submit/FAB1610H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

