

# Product Datasheet

## ICOS Antibody (2317C) [DyLight 550] FAB1682L

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB1682L](http://www.novusbio.com/FAB1682L)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB1682L](http://www.novusbio.com/reviews/destination/FAB1682L)



**FAB1682L**

ICOS Antibody (2317C) [DyLight 550]

| Product Information         |                                                                                              |
|-----------------------------|----------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                       |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.  |
| Storage                     | Store at 4C in the dark.                                                                     |
| Clonality                   | Monoclonal                                                                                   |
| Clone                       | 2317C                                                                                        |
| Preservative                | 0.05% Sodium Azide                                                                           |
| Isotype                     | IgG                                                                                          |
| Conjugate                   | DyLight 550                                                                                  |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                   |
| Buffer                      | 50mM Sodium Borate                                                                           |
| Product Description         |                                                                                              |
| Host                        | Rabbit                                                                                       |
| Gene ID                     | 29851                                                                                        |
| Species                     | Mouse                                                                                        |
| Specificity/Sensitivity     | Detects mouse ICOS in direct ELISAs.                                                         |
| Immunogen                   | Chinese Hamster Ovary cell line, CHO-derived mouse ICOS<br>Met1-Leu142<br>Accession # Q9WVS0 |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.            |
| Product Application Details |                                                                                              |
| Applications                | Flow Cytometry, Blockade of Receptor-ligand Interaction, CyTOF-ready                         |
| Recommended Dilutions       | Flow Cytometry, Blockade of Receptor-ligand Interaction, CyTOF-ready                         |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                       |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to FAB1682L**

---

|               |                                          |
|---------------|------------------------------------------|
| NBP2-24891R   | Rabbit IgG Isotype Control [DyLight 550] |
| NBP3-21257PEP | ICOS Recombinant Protein Antigen         |
| 210-TA-005    | TNF-alpha [Unconjugated]                 |
| DY169         | ICOS [Biotin]                            |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB1682L](http://www.novusbio.com/reviews/submit/FAB1682L)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

