

# Product Datasheet

## Kallikrein 8/Neuropsin Antibody (369910) [DyLight 594] FAB17191M

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB17191M](http://www.novusbio.com/FAB17191M)

Updated 5/28/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB17191M](http://www.novusbio.com/reviews/destination/FAB17191M)



**FAB17191M**

Kallikrein 8/Neurotrophin Antibody (369910) [DyLight 594]

| Product Information         |                                                                                                                                                                                                                                                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                          |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                     |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                        |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                      |
| Clone                       | 369910                                                                                                                                                                                                                                                                                          |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                              |
| Isotype                     | IgG1                                                                                                                                                                                                                                                                                            |
| Conjugate                   | DyLight 594                                                                                                                                                                                                                                                                                     |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                      |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                              |
| Product Description         |                                                                                                                                                                                                                                                                                                 |
| Host                        | Rat                                                                                                                                                                                                                                                                                             |
| Gene ID                     | 11202                                                                                                                                                                                                                                                                                           |
| Gene Symbol                 | KLK8                                                                                                                                                                                                                                                                                            |
| Species                     | Human                                                                                                                                                                                                                                                                                           |
| Specificity/Sensitivity     | Detects human Kallikrein 8/Neurotrophin in direct ELISAs. In direct ELISAs, no cross-reactivity with recombinant human (rh) UPA□1, rhKallikrein 1, 2, 3, 4, 5, 7, 9, 10, 11, 12, 13, 14, 15, rhKallikrein B1, rhHGFA, rhFactor VII, rhFactor X, rhFactor XI, rhProc, or rhThrombin is observed. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant human Kallikrein 8/Neurotrophin<br>Gln29-Gly260 (predicted)<br>Accession # O60259                                                                                                                                                               |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                               |
| Product Application Details |                                                                                                                                                                                                                                                                                                 |
| Applications                | Immunoprecipitation, Neutralization                                                                                                                                                                                                                                                             |
| Recommended Dilutions       | Immunoprecipitation, Neutralization                                                                                                                                                                                                                                                             |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                          |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to FAB17191M**

|                 |                                                        |
|-----------------|--------------------------------------------------------|
| NBP2-59574-50ug | Recombinant Mouse Kallikrein 8/Neurotensin His Protein |
| DKK300          | Kallikrein 3/PSA [HRP]                                 |
| 2025-SE-010     | Kallikrein 8/Neurotensin [Unconjugated]                |
| MAB208-100      | CXCL8/IL-8 Antibody (6217) [Unconjugated]              |

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB17191M](http://www.novusbio.com/reviews/submit/FAB17191M)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

