

Product Datasheet

MIS/AMH Antibody (213816) [DyLight 594] FAB1737M

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB1737M

Updated 5/28/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB1737M



FAB1737M

MIS/AMH Antibody (213816) [DyLight 594]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	213816
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	DyLight 594
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	50mM Sodium Borate
Product Description	
Host	Mouse
Gene ID	268
Gene Symbol	AMH
Species	Human
Specificity/Sensitivity	Detects human MIS/AMH in direct ELISAs and Western blots. Shows approximately 50-100% cross-reactivity with recombinant mouse (rm) MIS and recombinant rat MIS and no cross-reactivity with recombinant human (rh) ART, rmArtemin, rhCripto-1, recombinant drosophila Dpp, rhGDNF, rhLAP, rhLatent TGF- β 1, rhLefty-A (short and long isoforms), rmLefty-1 (short and long isoforms), rhNeurturin, rmNodal, rhTGF- α , recombinant amphibian TGF- β 5, rhTGF- β 1, rhTGF- β 1.2, rhTGF- β 2, rhTGF- β 3, or rmWIF-1.
Immunogen	<i>E. coli</i> -derived recombinant human MIS/AMH Ala453-Arg560 Accession # NP_000470
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB1737M

NBP1-97005DL594	Mouse IgG1 Isotype Control (MG1) [DyLight 594]
NBP2-68847PEP	MIS/AMH Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]
1737-MS-010	MIS/AMH [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB1737M

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

