

Product Datasheet

Reg3A Antibody (250802) [DyLight 550] FAB1745L

Unit Size: 0.5 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB1745L

Updated 7/14/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB1745L



FAB1745L

Reg3A Antibody (250802) [DyLight 550]

Product Information	
Unit Size	0.5 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	250802
Preservative	0.05% Sodium Azide
Isotype	IgG2a
Conjugate	DyLight 550
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	50mM Sodium Borate
Product Description	
Host	Mouse
Gene ID	5068
Gene Symbol	REG3A
Species	Rat
Specificity/Sensitivity	Detects rat Reg3A in direct ELISAs and Western blots. In these formats, shows 100% cross-reactivity with recombinant mouse (rm) Reg3A and no cross-reactivity with rmReg1, recombinant human (rh) Reg1B, recombinant rat Reg2, or rhReg4.
Immunogen	<i>E. coli</i> -derived recombinant rat Reg3A Glu26-Gln174 Accession # NP_742074
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB1745L

NBP2-56285PEP	Reg3A Recombinant Protein Antigen
DKK300	Kallikrein 3/PSA [HRP]
DY5940-05	Reg3A [Biotin]
NBP1-19773	Serpin E1/PAI-1 Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB1745L

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

