

Product Datasheet

IL-17RA/IL-17R Antibody (133617) [DyLight 550] FAB177L

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB177L

Updated 8/2/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB177L



FAB177L

IL-17RA/IL-17R Antibody (133617) [DyLight 550]

| Product Information | |
|---------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 133617 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 |
| Conjugate | DyLight 550 |
| Purity | Protein A or G purified from hybridoma culture supernatant |
| Buffer | 50mM Sodium Borate |

| Product Description | |
|-------------------------|---|
| Host | Mouse |
| Gene ID | 23765 |
| Gene Symbol | IL17RA |
| Species | Human |
| Specificity/Sensitivity | Detects human IL-17 RA/IL-17 R in direct ELISAs and Western blots. In direct ELISAs, no cross-reactivity with recombinant mouse (rm) IL-17 R, rmIL-17B R, or recombinant human IL-17B R is observed. In Western blots, no cross-reactivity with rmIL-17 R is observed |
| Immunogen | Mouse myeloma cell line NS0-derived recombinant human IL-17 RA/IL-17 R Leu33-Trp320 Accession # Q96F46 |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details | |
|-----------------------------|---|
| Applications | Western Blot, Flow Cytometry, CyTOF-ready, Knockout Validated, Neutralization |
| Recommended Dilutions | Western Blot, Flow Cytometry, Neutralization, CyTOF-ready, Knockout Validated |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB177L

| | |
|-------------|--|
| NBP1-97005R | Mouse IgG1 Isotype Control (MG1) [DyLight 550] |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| DY177 | IL-17RA/IL-17R [Biotin] |
| DY453-05 | CXCL1/GRO alpha/KC/CINC-1 [Biotin] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB177L

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

