

# Product Datasheet

## TREM2 Antibody (1057114) [PE] FAB18282P

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB18282P](http://www.novusbio.com/FAB18282P)

Updated 11/4/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB18282P](http://www.novusbio.com/reviews/destination/FAB18282P)



**FAB18282P**

TREM2 Antibody (1057114) [PE]

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C in the dark.                                                                    |
| Clonality                   | Monoclonal                                                                                  |
| Clone                       | 1057114                                                                                     |
| Preservative                | 0.05% Sodium Azide                                                                          |
| Isotype                     | IgG1                                                                                        |
| Conjugate                   | PE                                                                                          |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer                      | PBS                                                                                         |
| Product Description         |                                                                                             |
| Host                        | Mouse                                                                                       |
| Gene ID                     | 54209                                                                                       |
| Species                     | Human                                                                                       |
| Specificity/Sensitivity     | Detects Human Trem-2 in direct and capture ELISAs                                           |
| Immunogen                   | NS0-derived recombinant human TREM-2 protein<br>His19-Ser174                                |
| Product Application Details |                                                                                             |
| Applications                | Immunohistochemistry                                                                        |
| Recommended Dilutions       | Immunohistochemistry                                                                        |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                      |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to FAB18282P**

---

|              |                                       |
|--------------|---------------------------------------|
| NBP1-97005PE | Mouse IgG1 Isotype Control (MG1) [PE] |
| NBP2-22888   | Recombinant Human TREM2 His Protein   |
| 210-TA-005   | TNF-alpha [Unconjugated]              |
| 1729-T2-050  | TREM2 [Unconjugated]                  |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB18282P](http://www.novusbio.com/reviews/submit/FAB18282P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

