

# Product Datasheet

## CCR7 Antibody (150503) [DyLight 350] FAB197D

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB197D](http://www.novusbio.com/FAB197D)

Updated 7/30/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB197D](http://www.novusbio.com/reviews/destination/FAB197D)



**FAB197D**

CCR7 Antibody (150503) [DyLight 350]

| <b>Product Information</b> |                                                                                             |
|----------------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>           | 0.1 ml                                                                                      |
| <b>Concentration</b>       | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>             | Store at 4C in the dark.                                                                    |
| <b>Clonality</b>           | Monoclonal                                                                                  |
| <b>Clone</b>               | 150503                                                                                      |
| <b>Preservative</b>        | 0.05% Sodium Azide                                                                          |
| <b>Isotype</b>             | IgG2a                                                                                       |
| <b>Conjugate</b>           | DyLight 350                                                                                 |
| <b>Purity</b>              | Protein A or G purified from hybridoma culture supernatant                                  |
| <b>Buffer</b>              | 50mM Sodium Borate                                                                          |

| <b>Product Description</b>     |                                                                                   |
|--------------------------------|-----------------------------------------------------------------------------------|
| <b>Host</b>                    | Mouse                                                                             |
| <b>Gene ID</b>                 | 1236                                                                              |
| <b>Gene Symbol</b>             | CCR7                                                                              |
| <b>Species</b>                 | Human                                                                             |
| <b>Specificity/Sensitivity</b> | Detects human CCR7.                                                               |
| <b>Immunogen</b>               | Human CCR7 transfectants<br>Met1-Pro378<br>Accession # AAA58615                   |
| <b>Notes</b>                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| <b>Product Application Details</b> |                                                                                           |
|------------------------------------|-------------------------------------------------------------------------------------------|
| <b>Applications</b>                | Flow Cytometry, Immunohistochemistry, CyTOF-reported, Immunocytochemistry, Neutralization |
| <b>Recommended Dilutions</b>       | Flow Cytometry, Immunohistochemistry, Neutralization, Immunocytochemistry, CyTOF-reported |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.                    |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB197D**

---

|                |                                  |
|----------------|----------------------------------|
| NBP2-57375PEP  | CCR7 Recombinant Protein Antigen |
| 210-TA-005     | TNF-alpha [Unconjugated]         |
| NBP2-08171     | CCR7 Overexpression Lysate       |
| 6507-IL-010/CF | IL-4 [Unconjugated]              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB197D](http://www.novusbio.com/reviews/submit/FAB197D)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

