

# Product Datasheet

## TNF-alpha Antibody (6401) [DyLight 350] FAB2101D

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB2101D](http://www.novusbio.com/FAB2101D)

Updated 8/3/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB2101D](http://www.novusbio.com/reviews/destination/FAB2101D)



**FAB2101D**

TNF-alpha Antibody (6401) [DyLight 350]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	6401
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	DyLight 350
Purity	Protein A or G purified from ascites
Buffer	50mM Sodium Borate
Product Description	
Host	Mouse
Gene ID	7124
Gene Symbol	TNF
Species	Human
Specificity/Sensitivity	Detects human TNF- $\alpha$ in direct ELISAs and Western blots. In direct ELISAs, approximately 25-50% cross-reactivity with recombinant porcine TNF- $\alpha$ and recombinant rhesus macaque TNF- $\alpha$ is observed but no cross-reactivity with recombinant cotton rat TNF- $\alpha$ , recombinant rat TNF- $\alpha$ , recombinant human (rh) LT $\alpha$ 1/ $\beta$ 2, rhLT $\alpha$ 2/ $\beta$ 1, rhAPRIL, rhBAFF, rhEDA- $\Delta$ A2, recombinant mouse EDA, rhFas Ligand, rhLIGHT, rhOX40 Ligand, rhTRAIL, rhTRANCE, rhTWEAK, or rhVEGI is observed.
Immunogen	<i>E. coli</i> -derived recombinant human TNF- $\alpha$ Gly57-Leu233 (predicted) Accession # P01375
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Flow Cytometry, CyTOF-ready, Neutralization
Recommended Dilutions	Flow Cytometry, Neutralization, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB2101D**

---

NBP1-97005UV	Mouse IgG1 Isotype Control (MG1) [DyLight 350]
NBP2-35076-10ug	Recombinant Human TNF-alpha Protein
201-LB-005	IL-1 beta/IL-1F2 [Unconjugated]
210-TA-005	TNF-alpha [Unconjugated]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB2101D](http://www.novusbio.com/reviews/submit/FAB2101D)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

