

Product Datasheet

GITR Ligand/TNFSF18 Antibody (721926) [DyLight 680] FAB21772Y

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB21772Y

Updated 5/10/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB21772Y



FAB21772Y

GITR Ligand/TNFSF18 Antibody (721926) [DyLight 680]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	721926
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	DyLight 680
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	50mM Sodium Borate
Product Description	
Host	Rat
Gene ID	8995
Gene Symbol	TNFSF18
Species	Mouse
Specificity/Sensitivity	Detects mouse GITR Ligand/TNFSF18 in direct ELISAs and Western blots. In direct ELISAs, no cross-reactivity with recombinant mouse (rm) APRIL, rmBAFF, rmEDA, rmFas Ligand, rmLIGHT, rmOX40 Ligand, rmTNF- α , rmTRAIL, rmTRANCE, rmTWEAK, rm4-1BB Ligand, recombinant human (rh) EDA-A2, rhLTa1/b2, rhLTa2/b1, rhVEGI, or rhGITR Ligand is observed. In Western blots, no cross-reactivity with rm4-1BB, rmEDA, rmLIGHT, rmTRAIL, rhGITR Ligand, rhTNF- α , recombinant canine TNF- α , recombinant guinea pig TNF- α , or recombinant porcine TNF- α is observed.
Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse GITR Ligand/TNFSF18 Thr47-Ser173 Accession # NP_899247
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB21772Y

NBP2-26579	Recombinant Mouse GITR Ligand/TNFSF18 Protein
210-TA-005	TNF-alpha [Unconjugated]
2177-GL-025	GITR Ligand/TNFSF18
285-IF-100	IFN-gamma [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB21772Y

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

