

Product Datasheet

Histone H2AX [p Ser139] Antibody (2207D) [Allophycocyanin] FAB2288A

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB2288A

Updated 10/7/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB2288A



FAB2288A

Histone H2AX [p Ser139] Antibody (2207D) [Allophycocyanin]

| Product Information | |
|---------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 2207D |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG |
| Conjugate | Allophycocyanin |
| Purity | Protein A or G purified from cell culture supernatant |
| Buffer | PBS |

| Product Description | |
|-------------------------|---|
| Host | Rabbit |
| Gene ID | 3014 |
| Species | Human |
| Specificity/Sensitivity | Detects human Histone H2AX when phosphorylated at S139 in Western blots. |
| Immunogen | Phospho-peptide containing human Histone H2AX S139 site Accession # P16104 |

| Product Application Details | |
|-----------------------------|--|
| Applications | Western Blot, CyTOF-ready, Intracellular Staining by Flow Cytometry |
| Recommended Dilutions | Western Blot, Intracellular Staining by Flow Cytometry, CyTOF-ready |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB2288A

| | |
|------------|---|
| IC1051A | Rabbit IgG Isotype Control (60024B) [Allophycocyanin] |
| NBP1-42464 | Recombinant Human Histone H2AX His Protein |
| DYC2288-2 | Histone H2AX [p Ser139] [Biotin] |
| NB100-304 | 53BP1 Antibody - BSA Free |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB2288A

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

