

Product Datasheet

CD30/TNFRSF8 Antibody (81337) [HRP] FAB229H

Unit Size: 0.5 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB229H

Updated 7/30/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB229H



FAB229H**CD30/TNFRSF8 Antibody (81337) [HRP]**

| Product Information | |
|----------------------------|---|
| Unit Size | 0.5 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 81337 |
| Preservative | No Preservative |
| Isotype | IgG2b |
| Conjugate | HRP |
| Purity | Protein A or G purified from hybridoma culture supernatant |
| Buffer | PBS |

| Product Description | |
|--------------------------------|--|
| Host | Mouse |
| Gene ID | 943 |
| Gene Symbol | TNFRSF8 |
| Species | Human |
| Reactivity Notes | 0 |
| Specificity/Sensitivity | Detects human CD30/TNFRSF8 in direct ELISAs and Western blots. In direct ELISAs, less than 1% cross-reactivity with recombinant mouse CD30, recombinant human (rh) CD27, and rhCD40 is observed. |
| Immunogen | <i>S. frugiperda</i> insect ovarian cell line Sf 21-derived recombinant human CD30/TNFRSF8 Phe19-Lys379 Accession # P28908 |

| Product Application Details | |
|------------------------------------|--|
| Applications | Western Blot, Flow Cytometry, Agonist Activity, CyTOF-ready |
| Recommended Dilutions | Western Blot, Flow Cytometry, Agonist Activity, CyTOF-ready |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB229H

| | |
|-------------|--|
| NBP2-27231H | Mouse IgG2b Isotype Control (MPC-11) [HRP] |
| NBP2-22660 | Recombinant Human CD30/TNFRSF8 His Protein |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| 813-CD-100 | CD30/TNFRSF8 [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB229H

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

