

# Product Datasheet

## **TIM-3 Antibody (344801) [DyLight 488] FAB23651K**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB23651K](http://www.novusbio.com/FAB23651K)

Updated 5/28/2025 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB23651K](http://www.novusbio.com/reviews/destination/FAB23651K)



**FAB23651K**

TIM-3 Antibody (344801) [DyLight 488]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 344801                                                                                      |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG1                                                                                        |
| Conjugate           | DyLight 488                                                                                 |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description     |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| Host                    | Rat                                                                                              |
| Gene ID                 | 84868                                                                                            |
| Gene Symbol             | HAVCR2                                                                                           |
| Species                 | Human                                                                                            |
| Specificity/Sensitivity | Detects recombinant human TIM-3 in direct ELISAs and Western blots.                              |
| Immunogen               | Mouse myeloma cell line NS0-derived recombinant human TIM-3<br>Met1-Arg200<br>Accession # Q8TDQ0 |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Agonist Activity                                                       |
| Recommended Dilutions       | Agonist Activity                                                       |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to FAB23651K**

---

|                 |                                          |
|-----------------|------------------------------------------|
| NBP3-07123-50ug | Recombinant Mouse TIM-3 hlgG-His Protein |
| 210-TA-005      | TNF-alpha [Unconjugated]                 |
| 2365-TM-050     | TIM-3 [Unconjugated]                     |
| 6507-IL-010/CF  | IL-4 [Unconjugated]                      |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB23651K](http://www.novusbio.com/reviews/submit/FAB23651K)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

