

# Product Datasheet

## **DcTRAILR1/TNFRSF23 Antibody (270028) [Janelia Fluor® 669]**

**FAB23781JF669**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB23781JF669](http://www.novusbio.com/FAB23781JF669)

Updated 5/8/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB23781JF669](http://www.novusbio.com/reviews/destination/FAB23781JF669)



**FAB23781JF669**

DcTRAILR1/TNFRSF23 Antibody (270028) [Janelia Fluor® 669]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                          |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                              |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                                                 |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                                                                                               |
| Clone                       | 270028                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                                                                                       |
| Isotype                     | IgG2a                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Conjugate                   | Janelia Fluor 669                                                                                                                                                                                                                                                                                                                                                                                                        |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                                                                                                                                               |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                                                                                                       |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Host                        | Rat                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Gene ID                     | 8718                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Gene Symbol                 | Tnfrsf23                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Species                     | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Specificity/Sensitivity     | Detects mouse DcTRAIL R1/TNFRSF23 in direct ELISAs and Western blots. In Western blots, no cross-reactivity with recombinant human (rh) BAFF R, recombinant mouse (rm) BAFF R, rhCD27, rmCD27, rhCD30, rmCD30, rhCD40, rmCD40, rhDR3, rmDR3, rhDR6, rhEDAR, rmEDAR, rhFas, rmFas, recombinant rat Fas, recombinant feline Fas, rhGITR, rmGITR, rhLymphotoxin $\beta$ R, rhNGF R, rmNGF R, rhRELt, or rhTROY is observed. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant mouse DcTRAIL R1/TNFRSF23<br>Met30-Ser156<br>Accession # Q9ER63                                                                                                                                                                                                                                                                                                          |
| Notes                       | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.                                                                                                                                                                                                                                                                                                                                    |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Applications                | Western Blot                                                                                                                                                                                                                                                                                                                                                                                                             |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                                                                                                                                                                                                                                             |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                                                   |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to FAB23781JF669**

---

|             |                                   |
|-------------|-----------------------------------|
| 210-TA-005  | TNF-alpha [Unconjugated]          |
| 8479-DT-025 | DcTRAILR1/TNFRSF23 [Unconjugated] |
| DRT100      | TNF RI/TNFRSF1A [HRP]             |
| NBP1-19371  | CD4 Antibody - BSA Free           |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB23781JF669](http://www.novusbio.com/reviews/submit/FAB23781JF669)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

