

Product Datasheet

DcTRAILR1/TNFRSF23 Antibody (270012) [DyLight 488] FAB2378K

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB2378K

Updated 7/13/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB2378K



FAB2378K

DcTRAILR1/TNFRSF23 Antibody (270012) [DyLight 488]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	270012
Preservative	0.05% Sodium Azide
Isotype	IgG2a
Conjugate	DyLight 488
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	50mM Sodium Borate
Product Description	
Host	Rat
Gene ID	8718
Gene Symbol	Tnfrsf23
Species	Mouse
Specificity/Sensitivity	Detects mouse DcTRAIL R1/TNFRSF23 in direct ELISAs and Western blots. In Western blots, less than 10% cross-reactivity with recombinant mouse (rm) 4-1BB Ligand, recombinant human (rh) BAFF R, rhDcTRAIL R3, rmFas, rmOPG, and rhTRAIL R4 is observed and no cross-reactivity with rmCD27, rmCD30, rmDR3, rhDR6, rmEDAR, recombinant feline Fas, rhFas, recombinant rat Fas, rmGITR, rmHVEM, rmLymphotoxin beta R, rmNGF R, rmOX40, rmTNF RI, rmTNF RII, rhTRAIL R1, or rhTROY is observed.
Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse DcTRAIL R1/TNFRSF23 Met30-Ser156 Accession # Q9ER63
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot, Immunocytochemistry
Recommended Dilutions	Western Blot, Immunocytochemistry
Application Notes	Optimal dilution of this antibody should be experimentally determined.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB2378K

210-TA-005	TNF-alpha [Unconjugated]
8479-DT-025	DcTRAILR1/TNFRSF23 [Unconjugated]
DRT100	TNF RI/TNFRSF1A [HRP]
NBP1-19371	CD4 Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB2378K

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

