

Product Datasheet

NCAM-1/CD56 Antibody (301040) [HRP] FAB2408H

Unit Size: 0.5 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB2408H

Updated 7/29/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB2408H



FAB2408H**NCAM-1/CD56 Antibody (301040) [HRP]**

| Product Information | |
|----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size | 0.5 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 301040 |
| Preservative | No Preservative |
| Isotype | IgG2b |
| Conjugate | HRP |
| Purity | Protein A or G purified from hybridoma culture supernatant |
| Buffer | PBS |

| Product Description | |
|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host | Mouse |
| Gene ID | 4684 |
| Gene Symbol | NCAM1 |
| Species | Human, Primate |
| Reactivity Notes | 0 |
| Specificity/Sensitivity | Detects human NCAM \square 1/CD56 in ELISAs and Western blots. In direct ELISAs and Western blots, no cross-reactivity with recombinant human (rh) ALCAM, rhBCAM, rhEPCAM, rhMCAM, rhNCAM-1-L1, recombinant mouse (rm) MAdCAM-1, rmNCAM-1-L1, or rmOCAM is observed. |
| Immunogen | Mouse myeloma cell line NS0-derived recombinant human NCAM \square 1/CD56 |

| Product Application Details | |
|------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications | Western Blot, Flow Cytometry, CyTOF-ready, ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), ELISA Standard (Matched Pair) |
| Recommended Dilutions | Western Blot, Flow Cytometry, ELISA Standard (Matched Pair), ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), CyTOF-ready |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB2408H

| | |
|-----------------|--------------------------------------------|
| NBP2-27231H | Mouse IgG2b Isotype Control (MPC-11) [HRP] |
| NBP3-07122-50ug | Recombinant Human NCAM-1/CD56 His Protein |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| DY2408 | NCAM-1/CD56 [Biotin] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB2408H

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

