

Product Datasheet

TGF-beta 3 Antibody (20724) [DyLight 680] FAB243Y

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB243Y

Updated 7/13/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB243Y



FAB243Y

TGF-beta 3 Antibody (20724) [DyLight 680]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	20724
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	DyLight 680
Purity	Protein A or G purified from ascites
Buffer	50mM Sodium Borate
Product Description	
Host	Mouse
Gene ID	7043
Gene Symbol	TGFB3
Species	Multi-Species
Specificity/Sensitivity	Detects TGF- β 3 from multiple species in direct ELISAs and Western blots. In Western blots, less than 25% cross-reactivity with recombinant human (rh) TGF- β 1.2 and rhTGF- β 2 is observed, and less than 2% cross-reactivity with recombinant amphibian TGF- β 5 and recombinant human TGF- β 1 is observed. Neutralizes the biological activity of TGF- β 3 but not TGF- β 1, TGF- β 2, or TGF- β 5.
Immunogen	Spodoptera frugiperda, Sf 21 (baculovirus) derived recombinant human TGF- β 3 Ala301-Ser412 (Tyr340Phe) Accession # P10600
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot, Immunohistochemistry, Neutralization
Recommended Dilutions	Western Blot, Immunohistochemistry, Neutralization
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB243Y

NBP1-97005FR	Mouse IgG1 Isotype Control (MG1) [DyLight 680]
NBP2-52029-0.1mg	Recombinant Human TGF-beta 3 His Protein
7754-BH-005/CF	TGF-beta 1 [Unconjugated]
243-B3-002	TGF-beta 3 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB243Y

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

