

# Product Datasheet

## Osteoactivin/GPNMB Antibody (303802) [HRP] FAB2550H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB2550H](http://www.novusbio.com/FAB2550H)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB2550H](http://www.novusbio.com/reviews/destination/FAB2550H)



**FAB2550H**

Osteoactivin/GPNMB Antibody (303802) [HRP]

| Product Information         |                                                                                                                                                                                     |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                              |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                         |
| Storage                     | Store at 4C in the dark.                                                                                                                                                            |
| Clonality                   | Monoclonal                                                                                                                                                                          |
| Clone                       | 303802                                                                                                                                                                              |
| Preservative                | No Preservative                                                                                                                                                                     |
| Isotype                     | IgG2a                                                                                                                                                                               |
| Conjugate                   | HRP                                                                                                                                                                                 |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                          |
| Buffer                      | PBS                                                                                                                                                                                 |
| Product Description         |                                                                                                                                                                                     |
| Host                        | Mouse                                                                                                                                                                               |
| Gene ID                     | 10457                                                                                                                                                                               |
| Gene Symbol                 | GPNMB                                                                                                                                                                               |
| Species                     | Human                                                                                                                                                                               |
| Specificity/Sensitivity     | Detects human Osteoactivin/GPNMB in direct ELISAs and Western blots. In direct ELISAs and Western blots, no cross-reactivity with recombinant mouse Osteoactivin/GPNMB is observed. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant human Osteoactivin/GPNMB isoform 2<br>Ala22-Asn486<br>Accession # NP_002501                                                         |
| Product Application Details |                                                                                                                                                                                     |
| Applications                | Western Blot                                                                                                                                                                        |
| Recommended Dilutions       | Western Blot                                                                                                                                                                        |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                              |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to FAB2550H**

|            |                                                  |
|------------|--------------------------------------------------|
| NBP1-97075 | Mouse IgG2a Isotype Control (M2A) [HRP]          |
| NBP2-23082 | Recombinant Human Osteoactivin/GPNMB His Protein |
| DY2330     | Osteoactivin/GPNMB [Biotin]                      |
| AF808      | Osteopontin/OPN Antibody [Unconjugated]          |

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB2550H](http://www.novusbio.com/reviews/submit/FAB2550H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

