

# Product Datasheet

## HIF-2 alpha/EPAS1 Antibody (2597C) [PerCP] FAB29971C

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB29971C](http://www.novusbio.com/FAB29971C)

Updated 5/28/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB29971C](http://www.novusbio.com/reviews/destination/FAB29971C)



**FAB29971C**

HIF-2 alpha/EPAS1 Antibody (2597C) [PerCP]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 2597C                                                                                       |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG                                                                                         |
| Conjugate           | PerCP                                                                                       |
| Purity              | Protein A or G purified from cell culture supernatant                                       |
| Buffer              | PBS                                                                                         |

| Product Description     |                                                                                                          |
|-------------------------|----------------------------------------------------------------------------------------------------------|
| Host                    | Rabbit                                                                                                   |
| Gene ID                 | 2034                                                                                                     |
| Species                 | Human                                                                                                    |
| Specificity/Sensitivity | Detects mouse HIF-2 $\alpha$ /EPAS1 in direct ELISAs.                                                    |
| Immunogen               | <i>E. coli</i> -derived recombinant mouse HIF-2 $\alpha$ /EPAS1<br>Ser542-Thr874<br>Accession # AAH57870 |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Immunohistochemistry                                                   |
| Recommended Dilutions       | Immunohistochemistry                                                   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB29971C**

---

|                   |                                                          |
|-------------------|----------------------------------------------------------|
| NBP2-24891PCP     | Rabbit IgG Isotype Control [PerCP]                       |
| H00002034-P03-2ug | Recombinant Human HIF-2 alpha/EPAS1 GST (N-Term) Protein |
| DVE00             | VEGF [HRP]                                               |
| NB100-902HIF2     | HIF-2 alpha/EPAS1 Antibody Pack                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB29971C](http://www.novusbio.com/reviews/submit/FAB29971C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

