

# Product Datasheet

## LAR/PTPRF Antibody (384727) [DyLight 488] FAB3004K

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB3004K](http://www.novusbio.com/FAB3004K)

Updated 5/28/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB3004K](http://www.novusbio.com/reviews/destination/FAB3004K)



**FAB3004K**

LAR/PTPRF Antibody (384727) [DyLight 488]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	384727
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	DyLight 488
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	50mM Sodium Borate
Product Description	
Host	Rat
Gene ID	5792
Gene Symbol	PTPRF
Species	Human
Specificity/Sensitivity	Detects endogenous human LAR in Western blots. In Western blots, this antibody detects both full-length protein and the cleaved LAR extracellular domain (ECD). It does not cross-react with recombinant human PTPRM, PTPRK, DEP1, or PTPRγ.
Immunogen	Mouse myeloma cell line NS0-derived recombinant human LAR Ala27-Glu1251 Accession # NP_569707
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB3004K**

---

NBP1-81325PEP	LAR/PTPRF Recombinant Protein Antigen
291-G1-200	IGF-I/IGF-1 [Unconjugated]
9377-PF-050	LAR/PTPRF [Unconjugated]
AF114	CD45 Antibody [Unconjugated]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB3004K](http://www.novusbio.com/reviews/submit/FAB3004K)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

