

# Product Datasheet

## Keap1 Antibody (333116) [Allophycocyanin] FAB3024A

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB3024A](http://www.novusbio.com/FAB3024A)

Updated 5/28/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB3024A](http://www.novusbio.com/reviews/destination/FAB3024A)



**FAB3024A**

Keap1 Antibody (333116) [Allophycocyanin]

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C in the dark.                                                                    |
| Clonality                   | Monoclonal                                                                                  |
| Clone                       | 333116                                                                                      |
| Preservative                | 0.05% Sodium Azide                                                                          |
| Isotype                     | IgG2b                                                                                       |
| Conjugate                   | Allophycocyanin                                                                             |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer                      | PBS                                                                                         |
| Product Description         |                                                                                             |
| Host                        | Mouse                                                                                       |
| Gene ID                     | 9817                                                                                        |
| Gene Symbol                 | KEAP1                                                                                       |
| Species                     | Human, Mouse, Rat                                                                           |
| Specificity/Sensitivity     | Detects human, mouse and rat Keap1 in Western blots.                                        |
| Immunogen                   | <i>E. coli</i> -derived recombinant human Keap1<br>Ala90-Ile250<br>Accession # Q14145       |
| Product Application Details |                                                                                             |
| Applications                | Western Blot, Knockout Validated                                                            |
| Recommended Dilutions       | Western Blot, Knockout Validated                                                            |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                      |





### Novus Biologicals USA

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### Bio-Techne Canada

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### Bio-Techne Ltd

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### General Contact Information

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### Products Related to FAB3024A

|               |                                                        |
|---------------|--------------------------------------------------------|
| IC0041A       | Mouse IgG2b Isotype Control (133303) [Allophycocyanin] |
| NBP1-83106PEP | Keap1 Recombinant Protein Antigen                      |
| E3-310-050    | Keap1 [Unconjugated]                                   |
| H00008878-M01 | p62/SQSTM1 Antibody (2C11) - Azide and BSA Free        |

### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB3024A](http://www.novusbio.com/reviews/submit/FAB3024A)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

