

# Product Datasheet

## IL-17/IL-17A Antibody (1019404) [Allophycocyanin] FAB31711A

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB31711A](http://www.novusbio.com/FAB31711A)

Updated 5/28/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB31711A](http://www.novusbio.com/reviews/destination/FAB31711A)



**FAB31711A**

IL-17/IL-17A Antibody (1019404) [Allophycocyanin]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 1019404                                                                                     |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG2a                                                                                       |
| Conjugate           | Allophycocyanin                                                                             |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | PBS                                                                                         |

| Product Description     |                                                                                                                                                         |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Mouse                                                                                                                                                   |
| Gene ID                 | 3605                                                                                                                                                    |
| Species                 | Human                                                                                                                                                   |
| Specificity/Sensitivity | Detects human IL-17A and human IL-17A/IL-17F heterodimer in direct ELISAs.                                                                              |
| Immunogen               | Chinese Hamster Ovary cell line CHO-derived<br>Human IL17A<br>(Gly24Ala155)<br>Accession # Q16552<br>Human IL17F<br>(Arg31Gln163)<br>Accession # Q96PD4 |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Neutralization                                                         |
| Recommended Dilutions       | Neutralization                                                         |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB31711A**

---

|                |                                                       |
|----------------|-------------------------------------------------------|
| IC003A         | Mouse IgG2a Isotype Control (20102) [Allophycocyanin] |
| NBP2-35040-5ug | Recombinant Human IL-17/IL-17A Protein                |
| 210-TA-005     | TNF-alpha [Unconjugated]                              |
| DY421-05       | IL-17/IL-17A [Biotin]                                 |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB31711A](http://www.novusbio.com/reviews/submit/FAB31711A)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

