

Product Datasheet

CD31/PECAM-1 Antibody (377537) [PerCP] FAB33871C

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB33871C

Updated 9/11/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB33871C



FAB33871C**CD31/PECAM-1 Antibody (377537) [PerCP]**

| Product Information | |
|------------------------------------|--|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 377537 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 |
| Conjugate | PerCP |
| Purity | Protein A or G purified from hybridoma culture supernatant |
| Buffer | PBS |
| Product Description | |
| Host | Rat |
| Gene ID | 5175 |
| Gene Symbol | PECAM1 |
| Species | Porcine |
| Specificity/Sensitivity | Detects porcine CD31/PECAM \square 1 in direct ELISAs and Western blots. In direct ELISAs and Western blots, no cross-reactivity with recombinant human (rh) CD31, recombinant mouse (rm) CD31, rhVCAM-1, rhICAM-1, rhICAM-2, rhICAM-3, rmlCAM-5, or rmMADCAM-1 is observed. |
| Immunogen | Mouse myeloma cell line NS0-derived recombinant porcine CD31/PECAM \square 1 Gln28-Lys602 (predicted) Accession # Q95242 |
| Product Application Details | |
| Applications | Western Blot, Flow Cytometry, CyTOF-ready |
| Recommended Dilutions | Western Blot, Flow Cytometry, CyTOF-ready |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB33871C

| | |
|---------------|--|
| NBP2-54655PEP | CD31/PECAM-1 Recombinant Protein Antigen |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| 3628-PC-050 | CD31/PECAM-1 [Unconjugated] |
| DVE00 | VEGF [HRP] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB33871C

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

