

# Product Datasheet

## Cytokeratin 19 Antibody (1042608) [FITC] FAB35063F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB35063F](http://www.novusbio.com/FAB35063F)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB35063F](http://www.novusbio.com/reviews/destination/FAB35063F)



**FAB35063F**

Cytokeratin 19 Antibody (1042608) [FITC]

| Product Information |   |
|---------------------|---|
| Unit Size           | 0.1 ml  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Monoclonal  |
| Clone               | 1042608   |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG2b   |
| Conjugate           | FITC  |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | PBS   |

| Product Description     |   |
|-------------------------|---|
| Host                    | Mouse   |
| Gene ID                 | 3880  |
| Species                 | Human   |
| Specificity/Sensitivity | Detects human Cytokeratin 19 in direct ELISAs                                       |
| Immunogen               | <i>E. coli</i> -derived human Cytokeratin 19<br>Gln311-Lys370<br>Accession # P08727 |

| Product Application Details |   |
|-----------------------------|---|
| Applications                | Western Blot, Immunohistochemistry, Immunocytochemistry, Intracellular Staining by Flow Cytometry |
| Recommended Dilutions       | Western Blot, Immunohistochemistry, Intracellular Staining by Flow Cytometry, Immunocytochemistry |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                            |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB35063F**

---

|               |   |
|---------------|---|
| NBP2-27229    | Mouse IgG2b Isotype Control (MPC-11) [FITC] |
| NB100-687PEP  | Cytokeratin 19 Antibody Blocking Peptide    |
| 1129-ER-050   | ErbB2/Her2 [Unconjugated]                   |
| NBP1-78278PEP | Cytokeratin 19 Antibody Blocking Peptide    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB35063F](http://www.novusbio.com/reviews/submit/FAB35063F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

