

Product Datasheet

Integrin alpha 7 Antibody (334908) [PE/Atto594] FAB3518PEATT594

Unit Size: 0.1 ml

Store at 4C in the dark. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB3518PEATT594

Updated 7/30/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB3518PEATT594



FAB3518PEATT594

Integrin alpha 7 Antibody (334908) [PE/Atto594]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark. Do not freeze.
Clonality	Monoclonal
Clone	334908
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Conjugate	PE/Atto594
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	PBS
Product Description	
Host	Rat
Gene ID	3679
Gene Symbol	ITGA7
Species	Mouse
Specificity/Sensitivity	Detects mouse Integrin α 7 in direct ELISAs. Reactivity with other Integrin α 7 isoforms was not tested. In direct ELISAs, no cross-reactivity with recombinant mouse Integrin α 2, α 3, α 4, α 5, α 9, α E, α X, recombinant human Integrin α L, or α 11 is observed.
Immunogen	Chinese hamster ovary cell line CHO-derived recombinant mouse Integrin α 7 Phe34-Pro1036 Accession # NP_032424
Product Application Details	
Applications	Flow Cytometry
Recommended Dilutions	Flow Cytometry
Application Notes	Optimal dilution of this antibody should be experimentally determined. For optimal results using our Tandem dyes, please avoid prolonged exposure to light or extreme temperature fluctuations. These can lead to irreversible degradation or decoupling. When staining intracellular targets, specific attention to the fixation and permeabilization steps in your flow protocol may be required. Please contact our technical support team at technical@novusbio.com if you have any questions.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB3518PEATT594

H00003679-Q01-10ug	Recombinant Human Integrin alpha 7 GST (N-Term) Protein
210-TA-005	TNF-alpha [Unconjugated]
NBP1-86118PEP	Integrin alpha 7 Recombinant Protein Antigen
6507-IL-010/CF	IL-4 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB3518PEATT594

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

