

# Product Datasheet

## **CCL15/MIP-1 delta Antibody (88119) [DyLight 350] FAB3631D**

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB3631D](http://www.novusbio.com/FAB3631D)

Updated 7/13/2025 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB3631D](http://www.novusbio.com/reviews/destination/FAB3631D)



**FAB3631D**

CCL15/MIP-1 delta Antibody (88119) [DyLight 350]

| Product Information         |                                                                                                                                                                                                                                                                                                                                            |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.5 ml                                                                                                                                                                                                                                                                                                                                     |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                   |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                 |
| Clone                       | 88119                                                                                                                                                                                                                                                                                                                                      |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                         |
| Isotype                     | IgG1                                                                                                                                                                                                                                                                                                                                       |
| Conjugate                   | DyLight 350                                                                                                                                                                                                                                                                                                                                |
| Purity                      | Protein A or G purified from ascites                                                                                                                                                                                                                                                                                                       |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                         |
| Product Description         |                                                                                                                                                                                                                                                                                                                                            |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                      |
| Gene ID                     | 6359                                                                                                                                                                                                                                                                                                                                       |
| Gene Symbol                 | CCL15                                                                                                                                                                                                                                                                                                                                      |
| Species                     | Human                                                                                                                                                                                                                                                                                                                                      |
| Specificity/Sensitivity     | Detects human CCL15/MIP-1 $\delta$ in direct ELISAs and Western blots. In direct ELISAs, no cross-reactivity with recombinant human CCL1, 2, 3, 4, 5, 7, 8, 11, 13, 14, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26, recombinant mouse CCL1, 2, 3, 4, 5, 6, 7, 9, 11, 12, 17, 19, 20, 21, 22, 24, 25, or recombinant rat CCL20 is observed. |
| Immunogen                   | <i>E. coli</i> -derived recombinant human CCL15/MIP-1 $\delta$<br>Ser46-Ile113<br>Accession # Q16663                                                                                                                                                                                                                                       |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                          |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                            |
| Applications                | Western Blot, Neutralization                                                                                                                                                                                                                                                                                                               |
| Recommended Dilutions       | Western Blot, Neutralization                                                                                                                                                                                                                                                                                                               |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                     |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB3631D**

---

|                |                                                |
|----------------|------------------------------------------------|
| NBP1-97005UV   | Mouse IgG1 Isotype Control (MG1) [DyLight 350] |
| NBP2-35043-5ug | Recombinant Human CCL15/MIP-1 delta Protein    |
| DCP00          | CCL2/JE/MCP-1 [HRP]                            |
| DY363          | CCL15/MIP-1 delta [Biotin]                     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB3631D](http://www.novusbio.com/reviews/submit/FAB3631D)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

