

# Product Datasheet

## TIM-1/KIM-1/HAVCR Antibody (511227) [HRP] FAB3689H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB3689H](http://www.novusbio.com/FAB3689H)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB3689H](http://www.novusbio.com/reviews/destination/FAB3689H)



**FAB3689H**

TIM-1/KIM-1/HAVCR Antibody (511227) [HRP]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 511227                                                                                      |
| Preservative        | No Preservative                                                                             |
| Isotype             | IgG1                                                                                        |
| Conjugate           | HRP                                                                                         |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | PBS                                                                                         |

| Product Description     |                                                                                                                                                                                                                                                |
|-------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Mouse                                                                                                                                                                                                                                          |
| Gene ID                 | 26762                                                                                                                                                                                                                                          |
| Gene Symbol             | HAVCR1                                                                                                                                                                                                                                         |
| Species                 | Rat                                                                                                                                                                                                                                            |
| Specificity/Sensitivity | Detects rat TIM-1/KIM-1/HAVCR in ELISAs. In sandwich immunoassays, no cross-reactivity or interference with recombinant mouse (rm) TIM $\square$ 1, rmTIM-3, rmTIM-4, recombinant human (rh) TIM $\square$ 1, rhTIM-4, or rhNKp44 is observed. |
| Immunogen               | Mouse Myeloma cell line NS0-derived recombinant rat TIM-1/KIM-1/HAVCR Ser18-Val238<br>Accession # O54947                                                                                                                                       |

| Product Application Details |                                                                                                               |
|-----------------------------|---------------------------------------------------------------------------------------------------------------|
| Applications                | ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), ELISA Standard (Matched Pair) |
| Recommended Dilutions       | ELISA Standard (Matched Pair), ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair) |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                        |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to FAB3689H**

---

|                |                                             |
|----------------|---------------------------------------------|
| NBP1-97074     | Mouse IgG1 Isotype Control (MG1) [HRP]      |
| NB100-56421PEP | TIM-1/KIM-1/HAVCR Antibody Blocking Peptide |
| 210-TA-005     | TNF-alpha [Unconjugated]                    |
| DSKM100        | TIM-1/KIM-1/HAVCR [HRP]                     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB3689H](http://www.novusbio.com/reviews/submit/FAB3689H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

