

Product Datasheet

Aminopeptidase N/CD13 Antibody (986002) [FITC] FAB38152F

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB38152F

Updated 10/7/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB38152F



FAB38152F

Aminopeptidase N/CD13 Antibody (986002) [FITC]

Product Information

| | |
|----------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 986002 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2a |
| Conjugate | FITC |
| Purity | Protein A or G purified from hybridoma culture supernatant |
| Buffer | PBS |

Product Description

| | |
|--------------------------------|---|
| Host | Mouse |
| Gene ID | 290 |
| Species | Human |
| Specificity/Sensitivity | Detects human Aminopeptidase N/CD13. |
| Immunogen | Mouse myeloma cell line NS0-derived recombinant human Aminopeptidase N/CD13 Lys69-Lys967 Accession # AAA51719 |

Product Application Details

| | |
|------------------------------|--|
| Applications | Flow Cytometry, Immunohistochemistry, CyTOF-ready |
| Recommended Dilutions | Flow Cytometry, Immunohistochemistry, CyTOF-ready |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |

Images

Aminopeptidase N/CD13 Antibody (986002) [FITC] [FAB38152F] - Vial of FITC conjugated antibody. FITC is optimally excited at 498 nm by the Blue laser (488 nm) and has an emission maximum of 519 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB38152F

| | |
|---------------|---|
| NBP1-96848 | Mouse IgG2a Isotype Control (M2A) [FITC] |
| NBP2-33855PEP | Aminopeptidase N/CD13 Recombinant Protein Antigen |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| 3815-ZN-010 | Aminopeptidase N/CD13 [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB38152F

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

