

Product Datasheet

IGF-I R/IGF1R Antibody (33255) [DyLight 680] FAB391Y

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB391Y

Updated 7/30/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB391Y



FAB391Y

IGF-I R/IGF1R Antibody (33255) [DyLight 680]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	33255
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	DyLight 680
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	50mM Sodium Borate
Product Description	
Host	Mouse
Gene ID	3480
Gene Symbol	IGF1R
Species	Human, Mouse
Specificity/Sensitivity	Detects human IGF-I R/IGF1R in sandwich ELISAs and Western blots. Detects mouse IGF-I R/IGF1R in Immunohistochemistry. In sandwich immunoassays, less than 0.15% cross-reactivity or interference was observed with recombinant human (rh) IGF-I, rhIGF-1, rhIL-3 R α , rhIL-9 R, and rhTGF- β RII.
Immunogen	<i>S. frugiperda</i> insect ovarian cell line Sf 21-derived recombinant human IGF-I R/IGF1R Glu31-Asn932 Accession # P08069
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunohistochemistry, ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), ELISA Standard (Matched Pair), Immunocytochemistry, Neutralization
Recommended Dilutions	Western Blot, Flow Cytometry, Immunohistochemistry, ELISA Standard (Matched Pair), Neutralization, ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), Immunocytochemistry
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB391Y

NBP1-97005FR	Mouse IgG1 Isotype Control (MG1) [DyLight 680]
NBP1-77680PEP	IGF-I R/IGF1R Antibody Blocking Peptide
236-EG-200	EGF [Unconjugated]
391-GR-050	IGF-I R/IGF1R [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB391Y

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

