

Product Datasheet

TAPP1/PLEKHA1 Antibody (376409) [HRP] FAB4039H

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB4039H

Updated 5/8/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB4039H



FAB4039H

TAPP1/PLEKHA1 Antibody (376409) [HRP]

| Product Information | |
|-----------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 376409 |
| Preservative | No Preservative |
| Isotype | IgG2a |
| Conjugate | HRP |
| Purity | Protein A or G purified from hybridoma culture supernatant |
| Buffer | PBS |
| Product Description | |
| Host | Rat |
| Gene ID | 59338 |
| Gene Symbol | PLEKHA1 |
| Species | Human |
| Specificity/Sensitivity | Detects human TAPP1 in direct ELISAs and Western blots. |
| Immunogen | <i>E. coli</i> -derived recombinant human TAPP1 Pro2-Val404 Accession # Q9HB21 |
| Product Application Details | |
| Applications | Western Blot |
| Recommended Dilutions | Western Blot |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB4039H

| | |
|--------------------|---|
| H00059338-P01-10ug | Recombinant Human TAPP1/PLEKHA1 GST (N-Term) Protein |
| NBP1-86967PEP | TAPP1/PLEKHA1 Recombinant Protein Antigen |
| AF887 | AKT [p Ser473] Antibody [Unconjugated] - Pan Specific |
| NB200-540 | Complement C3 Antibody (11H9) - BSA Free |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB4039H

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

